

Probe Compatibility Chart

Probe	Order Code	Advanced Gage Bodies			Standard Gage Bodies		Other			
		PosiTector Advanced Gage Body (S/N 784,000+)	PosiTector Advanced Gage Body (S/N 730,000 - 783,999)	PosiTector Advanced Gage Body (S/N 700,000 - 729,999)	PosiTector Standard Gage Body w/ Color Display (S/N 795,816+)	PosiTector Standard Gage Body w/ Non-Color Display (S/N 700,000 - 795,815)	PosiTector SmartLink	PosiTector 6000 Gage Bodies (S/N 600,000 - 699,999)	PosiTector UTG Gage Bodies (S/N 600,000 - 699,999)	PosiTector DPM Gage Bodies (S/N 600,000 - 699,999)
PosiTector 6000 F Probe	PRBF-D	X	X	X	X	X	X	X		
PosiTector 6000 FS Probe	PRBFS-D	X	X	X	X	X	X	X		
PosiTector 6000 FXS Probe	PRBFXS-D	X	X	X	X	X	X*			
PosiTector 6000 FRS Probe	PRBFRS-D	X	X	X	X	X	X	X		
PosiTector 6000 FOS Probe	PRBFOS-D	X	X	X	X	X	X	X		
PosiTector 6000 F45S Probe	PRBF45S-D	X	X	X	X	X	X	X		
PosiTector 6000 F90S Probe	PRBF90S-D	X	X	X	X	X	X	X		
PosiTector 6000 FT Probe	PRBFT-D	X	X	X	X	X	X	X		
PosiTector 6000 FTS Probe	PRBFTS-D	X	X	X	X	X	X	X		
PosiTector 6000 FTRS Probe	PRBFTRS-D	X	X	X	X	X	X	X		
PosiTector 6000 FHXS Probe	PRBFHXS-D	X	X	X	X	X	X			
PosiTector 6000 FKS Probe	PRBFKS-D	X	X	X	X	X	X	X		
PosiTector 6000 FLS Probe	PRBFLS-D	X	X	X	X	X	X			
PosiTector 6000 N Probe	PRBN-D	X	X	X	X	X	X	X		
PosiTector 6000 NS Probe	PRBNS-D	X	X	X	X	X	X	X		
PosiTector 6000 NRS Probe	PRBNRS-D	X	X	X	X	X	X	X		
PosiTector 6000 NAS Probe	PRBNAS-D	X	X	X	X	X	X	X		
PosiTector 6000 NOS Probe	PRBNOS-D	X	X	X	X	X	X	X		
PosiTector 6000 N45S Probe	PRBN45S-D	X	X	X	X	X	X	X		
PosiTector 6000 N90S Probe	PRBN90S-D	X	X	X	X	X	X	X		
PosiTector 6000 NKS Probe	PRBNKS-D	X	X	X	X	X	X	X		
PosiTector 6000 FN Probe	PRBFN-D	X	X	X	X	X	X	X		
PosiTector 6000 FNS Probe	PRBFNS-D	X	X	X	X	X	X	X		
PosiTector 6000 FNDS Probe	PRBFNDS-D	X	X		X					
PosiTector 6000 FNRS Probe	PRBFNRS-D	X	X	X	X	X	X	X		
PosiTector 6000 FNFS Probe	PRBFNFS-D	X	X	X	X	X	X	X		
PosiTector 6000 FNGS Probe	PRBFNGS-D	X	X	X	X	X	X			
PosiTector 200 B Probe	PRB200B-C	X	X	X	X	X				
PosiTector 200 C Probe	PRB200C-C	X	X	X	X	X				
PosiTector 200 D Probe	PRB200D-C	X	X	X	X	X				
PosiTector SPG Probe	PRBSPG-B	X	X	X	X	X	X			
PosiTector SPGS Probe	PRBSPGS-B	X	X	X	X	X	X			
PosiTector SPGCS Probe	PRBSPGCS-B	X	X	X	X					
PosiTector SPGTS Probe	PRBSPGTS-B	X	X	X	X					
PosiTector RTR H Probe	PRBRTRH	X	X	X	X	X	X			
PosiTector RTR 3D Probe	PRBRTR3D	X			X					
PosiTector SHD D Probe	PRBSHDD	X	X		X					
PosiTector SHD A Probe	PRBSHDA	X	X		X					
PosiTector BHI Probe	PRBBHI	X	X		X					
PosiTector UTG C Probe	PRBUTGC-C	X	X	X	X	X				
PosiTector UTG CX Probe	PRBUTGCX-C	X	X	X	X	X				
PosiTector UTG M Probe	PRBUTGM-C	X	X	X	X	X				
PosiTector UTG STD	N/A								X	
PosiTector UTG ME	N/A								X	
PosiTector DPM Probe	PRBDPM-C	X	X	X	X	X	X			
PosiTector DPMIR Probe	PRBDPMIR-C	X	X	X	X					
PosiTector DPMS Probe	PRBDPMS-C	X	X	X	X	X	X			
PosiTector DPM Probe (2-Stem)	PRBDPM-B	X	X	X	X	X	X			X
PosiTector DPMS Probe (2-Stem)	PRBDPMS-B	X	X	X	X	X	X			X
PosiTector DPMA Probe	PRBDPMA-C	X	X		X					
PosiTector DPMD Probe	PRBDPMD-C	X	X		X					
PosiTector IRT Probe	PRBIRT	X	X		X					
PosiTector CMM IS Probe	CMMIS1PK	X								
PosiTector SST Probe	PRBSST	X	X		X					

* PosiTector 6000 FXS only compatible with SmartLink has been updated.